



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Wide-neck screw cap bottle 1000 ml, amber glass, DIN 68, without screw cap

€6.80 As low as

€5.69

**plus VAT &
Shipping**

Product Images



Description

Wide neck packaging bottles made of brown soda glass, round, according to DIN.
Accessories: Screw thread caps according to DIN (Order No. 2-1495 ff)

Additional Information

No.	DT-0086
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1484
EAN	4058072228620
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	1000 ml
Color bottle	brown
Height	180.5 mm
Material	Soda-lime glass (normal glass) 2
DM bottle	103 mm
DM Neck	68 mm
Thread	DIN 68
Neck size	far
Material bottle	Soda-lime glass (normal glass) 2

