



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Wide-neck screw cap bottle 50 ml,
clear glass, DIN 40, without screw cap**

€1.26 As low as
€1.05
**plus VAT &
Shipping**

Product Images



Description

Wide neck bottles made of clear soda glass, round, acc. to DIN.
Accessories: Screw thread caps acc. to DIN (order no. 2-1495 ff.)

Additional Information

No.	DT-0076
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1471
EAN	4058072228521
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	50 ml
Color bottle	transparent (clear)
Height	69.5 mm
Material	Soda-lime glass (normal glass) 2
DM bottle	44 mm
DM Neck	40 mm
Thread	DIN 40
Neck size	far
Material bottle	Soda-lime glass (normal glass) 2

