



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Wide-neck screw cap bottle 250 ml,  
clear glass, DIN 55, without screw cap**

**€2.16** As low as  
**€1.80**  
**plus VAT &  
Shipping**



## Description

---

Wide neck bottles made of clear soda glass, round, acc. to DIN.  
Accessories: Screw thread caps acc. to DIN (order no. 2-1495 ff.)

## Additional Information

---

No.	DT-0078
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1473
EAN	4058072228545
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	250 ml
Color bottle	transparent (clear)
Height	113 mm
Material	Soda-lime glass (normal glass) 2
Thread	DIN 55
Neck size	far
DM bottle	63 mm
DM Neck	55 mm
Material bottle	Soda-lime glass (normal glass) 2

