



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Wide-neck screw cap bottle 500 ml, clear glass, DIN 55, without screw cap

**€3.92** As low as  
**€3.28**  
**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Wide neck bottles made of clear soda glass, round, acc. to DIN.  
Accessories: Screw thread caps acc. to DIN (order no. 2-1495 ff.)

## Additional Information

---

|                        |                                  |
|------------------------|----------------------------------|
| No.                    | DT-0079                          |
| Manufacturer (Brand)   | neoLab                           |
| old neoLab article no. | 2-1474                           |
| EAN                    | 4058072228552                    |
| Transport temperature  | Room temperature                 |
| autoclavable at 121°C  | Yes                              |
| Color                  | transparent (clear)              |
| Material               | Soda-lime glass (normal glass) 2 |
| Volume                 | 500 ml                           |
| Color bottle           | transparent (clear)              |
| Height                 | 154 mm                           |
| Material               | Soda-lime glass (normal glass) 2 |
| DM bottle              | 84 mm                            |
| DM Neck                | 55 mm                            |
| Thread                 | DIN 55                           |
| Neck size              | far                              |
| Material bottle        | Soda-lime glass (normal glass) 2 |

