



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PET wide neck bottles clear, 250 ml, 10 pcs./pack

**€23.50
plus VAT &
Shipping**

Product Images



Description

Wide-neck bottles made of crystal clear thermoplastic polyester
Closures made of white PP.

Food safe, high chemical resistance. Largely diffusion resistant.
Temperature range: -40°C to +60°C

Additional Information

No.	3-2101
Manufacturer (Brand)	neoLab
EAN	4058072018962
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polyethylene terephthalate (PET)
Volume	250 ml
sterile	No
Color bottle	transparent (crystal clear)
Height	75 mm
Type closure/ lid	Screw cap
TBST MAX	60 °C
TBST MIN	-40 °C
Bottles properties	with closure/lid closure autoclavable
Material	Polyethylene terephthalate (PET)
Neck size	far
Closure / lid properties	closed
DM bottle	80 mm
Color closure/ lid	white
Material bottle	Polyethylene terephthalate (PET)
Material closure/ lid	Polypropylene (PP)

