



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PET wide neck bottles clear, 500 ml, 10 pcs./pack

€34.00
plus VAT &
Shipping

Product Images



Description

Wide-neck bottles made of crystal clear thermoplastic polyester
Closures made of white PP.

Food safe, high chemical resistance. Largely diffusion resistant.
Temperature range: -40°C to +60°C

Additional Information

No.	3-2102
Manufacturer (Brand)	neoLab
EAN	4058072018979
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polyethylene terephthalate (PET)
Volume	500 ml
sterile	No
Color bottle	transparent (crystal clear)
Height	110 mm
Type closure/ lid	Screw cap
TBST MAX	60 °C
TBST MIN	-40 °C
Bottles properties	with closure/lid closure autoclavable
Material	Polyethylene terephthalate (PET)
DM bottle	95 mm
Color closure/ lid	white
Neck size	far
Material bottle	Polyethylene terephthalate (PET)
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed

