



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® PET wide neck bottles clear, 3000 ml,  
10 pcs./pack**

**€125.50  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Wide-neck bottles made of crystal clear thermoplastic polyester  
Closures made of white PP.

Food safe, high chemical resistance. Largely diffusion resistant.  
Temperature range: -40°C to +60°C

## Additional Information

No.	3-2106
Manufacturer (Brand)	neoLab
EAN	4058072019013
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polyethylene terephthalate (PET)
Volume	3000 ml
sterile	No
Color bottle	transparent (crystal clear)
Height	210 mm
Type closure/ lid	Screw cap
TBST MAX	60 °C
TBST MIN	-40 °C
Bottles properties	with closure/lid closure autoclavable
Material	Polyethylene terephthalate (PET)
Neck size	far
Closure / lid properties	closed
DM bottle	150 mm
Color closure/ lid	white
Material bottle	Polyethylene terephthalate (PET)
Material closure/ lid	Polypropylene (PP)

