



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Yutility bottle, 500 ml, 4 pcs./pack

€50.64
plus VAT &
Shipping

Product Images



Description

Specially shaped recessed grips ensure easier and safer handling. The optimized thread enables faster opening and closing of the bottle and is compatible with all GL 45 closures and other accessories. Space-saving use due to the slim shape when used in autoclaves up to 121°C and laboratory refrigerators.

Labeling area marked on the bottles perfectly matches the self-adhesive DURAN® YOUTILITY labels (Order No. 2-3041). Optimized scale shows the nominal volume of the bottle, allowing faster identification of the total volume.

Made of DURAN® borosilicate glass 3.3. (glass type I). The bottle features very good chemical resistance and high temperature resistance.

The printed retrace code enables the bottle to be traced.

Each DURAN® YOUTILITY bottle is supplied as a complete system with a pouring ring (cyan, polypropylene) and a GL 45 screw cap (cyan, polypropylene). The ergonomically designed screw cap with optimized grooves and ribs enables simplified, faster opening and closing of the bottle, especially for users with small hands or when working with gloves.

The labeling area marked on the screw cap perfectly matches the self-adhesive DURAN® YOUTILITY labels.

Additional Information

No.	SD-0920
Manufacturer (Brand)	DURAN
old neolab article no.	2-3045
EAN	4058072346386
Transport temperature	Room temperature
Color	transparent (clear)
Volume	500 ml
sterile	No
Color bottle	transparent (clear)
Height	193 mm
Type closure/ lid	Screw cap
TBST MAX	150 °C
TBST MIN	-40 °C
Bottles properties	with high hydrolytic resistance with closure/lid with recessed grip with batch certificate with labeling field with pouring ring graduated cap autoclavable bottle autoclavable
Thread	GL 45
Neck size	narrow
Closure / lid properties	closed
DM bottle	78 mm
DM Neck	45 mm
Spare part	SD-2731 SD-2755
Color closure/ lid	blue
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)

