




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Two-neck dosing bottle, HDPE, 1000 ml

**€4.57**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **dosing bottles with the two-neck system** are ideally suited for dosing a small amount of liquid from a storage vessel into a second chamber. This allows, for example, electrodes to be calibrated directly in the vessel without affecting the stock solution.

The graduation in the dosing chamber is clearly visible for 10, 20, 30, 40, 50, 60 ml.

- made of HDPE
- suitable tamper-evident closure: 2-1478

## Additional Information

---

No.	2-1477
Manufacturer (Brand)	neoLab
EAN	4058072155766
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume	1000 ml
Wide	132 mm
Height	253 mm
Depth	64 mm
Type closure/ lid	Screw cap
Material	High-density polyethylene (HDPE)
Floor outside	flat
Floor inside	flat
Thread	DIN 28
Neck size	narrow
Closure / lid properties	closed
Color closure/ lid	white

