



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® aluminum bottle 625 ml with UN approval, 73 x 206 mm

**€10.00** As low as  
**€8.58**  
plus VAT &  
Shipping

### Product Images



## Description

---

Made of 99.5% pure aluminum, seamlessly manufactured according to DIN EN 573.

Aluminum bottle System 51 UN

Unrestrictedly approved for the transport of dangerous goods.

UN-approved, tested for leak tightness and marked with test code.

Silver matt pickled inside and outside.

PE male thread with smooth neck passage and PE screw cap.

Packing group I (X)

## Additional Information

No.	2-2262
Manufacturer (Brand)	neoLab
EAN	4058072060206
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Volume	625 ml
sterile	No
Color bottle	silver
Height	206 mm
Type closure/ lid	Screw cap
Bottles properties	with special thermal shock resistance with closure/lid with seal with batch certificate bottle autoclavable
Material	Aluminum
Neck size	narrow
Material gasket / septum	High-density polyethylene (HDPE)
Closure / lid properties	closed with seal
DM bottle	73 mm
DM Neck	42 mm
Color closure/ lid	silver
Material bottle	Aluminum
Material closure/ lid	High-density polyethylene (HDPE)

