



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® aluminum bottle 1250 ml with UN approval, 88 x 256 mm

€15.00 As low as
€12.40
plus VAT &
Shipping

Product Images



Description

Made of 99.5% pure aluminum, seamlessly manufactured according to DIN EN 573.

Aluminum bottle System 51 UN

Unrestrictedly approved for the transport of dangerous goods.

UN-approved, tested for leak tightness and marked with test code.

Silver matt pickled inside and outside.

PE male thread with smooth neck passage and PE screw cap, with PE foam seal.

Packing group I (X)

Additional Information

No.	2-2263
Manufacturer (Brand)	neoLab
EAN	4058072060213
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Volume	1250 ml
sterile	No
Color bottle	silver
Height	256 mm
Type closure/ lid	Screw cap
Bottles properties	with special thermal shock resistance with closure/lid with seal with batch certificate bottle autoclavable
Material	Aluminum
Neck size	narrow
Material gasket / septum	High-density polyethylene (HDPE)
Closure / lid properties	closed with seal
DM bottle	88 mm
DM Neck	42 mm
Color closure/ lid	silver
Material bottle	Aluminum
Material closure/ lid	High-density polyethylene (HDPE)

