



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® aluminum bottle 2500 ml with UN approval, 120 x 276 mm

€26.70 As low
as
€22.90
plus VAT &
Shipping

Product Images



Description

Made of 99.5% pure aluminum, seamlessly manufactured according to DIN EN 573.

Aluminum bottle System 35 UN

Unrestrictedly approved for the transport of dangerous goods.

UN-approved, tested for leak tightness and marked with test code.

Silver matt pickled inside and outside.

PE external thread with smooth neck passage and PE screw cap, self-sealing, with PE foam seal.

Packing group I (X)

Additional Information

No.	2-2264
Manufacturer (Brand)	neoLab
EAN	4058072060220
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Volume	2500 ml
sterile	No
Color bottle	silver
Height	276 mm
Type closure/ lid	Screw cap
Bottles properties	with special thermal shock resistance with closure/lid with seal with batch certificate bottle autoclavable
Material	Aluminum
Neck size	narrow
Material gasket / septum	High-density polyethylene (HDPE)
Closure / lid properties	closed with seal
DM bottle	120 mm
DM Neck	42 mm
Color closure/ lid	silver
Material bottle	Aluminum
Material closure/ lid	High-density polyethylene (HDPE)

