



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® aluminum bottle 1250 ml, UN approval, neck 26 mm

€28.30
plus VAT &
Shipping

Product Images



Description

Bottles made of pure aluminum (Al 99.5% according to DIN 1712) with UN approval, seamlessly manufactured.
Silver matt pickled on the outside and inside.
Aluminum screw stopper with rubber ring.
Each bottle is individually tested for leak tightness.
High light protection function, low weight and high durability.
Packaging group I (X), type-tested, unsupported.

Additional Information

| | |
|--------------------------|---|
| No. | 7-3017 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072038472 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Aluminum |
| Volume | 1250 ml |
| sterile | No |
| Color bottle | silver |
| Height | 200 mm |
| Type closure/ lid | Screw cap |
| Bottles properties | with special thermal shock resistance with cap/lid cap autoclavable bottle autoclavable |
| Material | Aluminum |
| DM bottle | 100 mm |
| DM Neck | 26 mm |
| Color closure/ lid | silver |
| Neck size | far |
| Material gasket / septum | High-density polyethylene (HDPE) |
| Material bottle | Aluminum |
| Material closure/ lid | Aluminum |
| Closure / lid properties | closed with seal |

