



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® aluminum bottle 125 ml, UN approval, neck 16 mm

**€13.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Bottles made of pure aluminum (Al 99.5% according to DIN 1712) with UN approval, seamlessly manufactured.  
Silver matt pickled on the outside and inside.  
Aluminum screw stopper with rubber ring.  
Each bottle is individually tested for leak tightness.  
High light protection function, low weight and high durability.  
Packaging group I (X), type-tested, unsupported.

## Additional Information

No.	7-3018
Manufacturer (Brand)	neoLab
EAN	4058072038489
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Volume	125 ml
sterile	No
Color bottle	silver
Height	110 mm
Type closure/ lid	Screw cap
Bottles properties	with special thermal shock resistance with cap/lid cap autoclavable bottle autoclavable
Material	Aluminum
DM bottle	50 mm
DM Neck	16 mm
Color closure/ lid	silver
Neck size	narrow
Material gasket / septum	High-density polyethylene (HDPE)
Material bottle	Aluminum
Material closure/ lid	Aluminum
Closure / lid properties	closed with seal

