



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® feeding bottle with glass olive, Erlenmeyer shape, 1000 ml

€76.96 As low
as
€59.20
plus VAT &
Shipping

Product Images



Description

DURAN® Vacuum flask, with glass olive, Erlenmeyer shape, 1000 ml A vacuum flask, which is suitable for vacuum applications due to its thick wall thickness, is used as a collecting vessel for the filtrate during vacuum filtration. The vacuum is applied via a side connection (tube, glass olive or plastic olive). The large neck opening is used to accommodate the rubber sleeve as well as the filter chute. Suction flasks in small volumes (< 2000 ml) are available in an Erlenmeyer shape, suction flasks with a volume > 2000 ml have a bottle shape. Special features and advantages: The design of suction flasks with a tube has not only simplified and facilitated work in preparative and analytical laboratories, but at the same time reduced the risk of accidents Vacuum-proof due to the solid wall thickness Not in accordance with the Equipment and Product Safety Act Burnt-in and therefore permanently durable white printed image Erlenmeyer shape Application: Suction of liquids or vacuum filtration, Filtration at reduced pressure, Filtration, Vacuum filtration, Safe working under vacuum

Additional Information

No.	SD-0358
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112120154
EAN	4058072177270
Manufacturer EAN	4032051006041
Gross weight	0.848
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.77000 kg
Material	Borosilicate glass 3.3
Volume	1000 ml
Country of origin	Germany
sterile	No
Region of origin	Rhineland-Palatinate
Packaging volume	0.04222200 cc
Customs tariff number	70172000

