




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® feeding bottle with KECK mounting set, bottle shape, 10000 ml

€364.95
plus VAT &
Shipping

Product Images



Description

DURAN® Vacuum flask with KECK™ mounting set, bottle shape, 10000 ml A vacuum flask, which is suitable for vacuum applications due to its thick wall thickness, is used as a collection vessel for the filtrate during vacuum filtration. The vacuum is applied via a side connection (tube, glass olive or plastic olive). The large neck opening is used to accommodate the rubber sleeve as well as the filter chute. Suction flasks in small volumes (< 2000 ml) are available in an Erlenmeyer shape, suction flasks with a volume > 2000 ml have a bottle shape. Special features and advantages: Vacuum-proof due to solid wall thickness These suction flasks comply with the regulations of the Equipment and Product Safety Act The plastic olives can be replaced Burned-in and therefore permanently durable white printed image Bottle shape Application: Suction of liquids or vacuum filtration, filtration at reduced pressure, filtration, vacuum filtration, safe working under vacuum

Additional Information

No.	SD-0352
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112119486
EAN	4058072341824
Manufacturer EAN	4032051005914
Gross weight	5.575
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	4.76500 kg
Material	Borosilicate glass 3.3
Volume	10000 ml
Country of origin	Germany
sterile	No
Region of origin	Baden-Württemberg
Packaging volume	0.02629100 cc
Customs tariff number	70172000

