



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine replacement aspirate bottle 4 l with screw cap and olives

€78.04
plus VAT &
Shipping

Product Images



Description

Replacement suction bottle with screw cap and olives

Plastic suction bottle, vacuum-proof.

Order No. 2-9338: 4 liter capacity

Accessories: Quick coupling (Order No. 2-1868 and 2-1871)

Accessories for suction device (Order No. 2-9335 and 2-9336)

Additional Information

No.	PK-29338
Manufacturer (Brand)	neoLabLine
EAN	4058072010638
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
sterile	No
Type closure/ lid	Screw cap with olives
Bottles properties	with closure/lid
for DM hose inside 1	6.0 mm
Material collection bottle	Polypropylene (PP)
Volume collection bottle	4 l

