



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® feeding bottle, 500 ml Erlenmeyer shape, with mounting set

€53.60 As low
as
€41.23
plus VAT &
Shipping

Product Images



Description

Erlenmeyer-shaped suction flasks for filtration applications with magnetic filter holders (8-3000 ff) or filter chutes. The plastic olives are interchangeable. The suction flasks are vacuum-proof and comply with the regulations of the Equipment and Product Safety Act, DIN 12 476, ISO 6556.

Additional Information

No.	SD-0367
Manufacturer (Brand)	DURAN
old neoLab article no.	8-3012
EAN	4058072176891
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	105 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
Form	trapezoidal
Height	185 mm
Spare part	SD-2793
Type connection	Screw connection
Type glass jar	Feed bottle Erlenmeyer form

