



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® feeding bottle, 2000 ml Erlenmeyer shape, with mounting set

€127.04
plus VAT &
Shipping

Product Images



Description

Erlenmeyer-shaped suction flasks for filtration applications with magnetic filter holders (8-3000 ff) or filter chutes. The plastic olives are interchangeable. The suction flasks are vacuum-proof and comply with the regulations of the Equipment and Product Safety Act, DIN 12 476, ISO 6556.

Additional Information

No.	SD-0369
Manufacturer (Brand)	DURAN
old neoLab article no.	8-3014
EAN	4058072341947
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	166 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	2000 ml
Form	trapezoidal
Height	255 mm
Spare part	SD-2793
Type connection	Screw connection
Type glass jar	Feed bottle Erlenmeyer form

