



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Unitary™ LDPE Squirt Bottles 250 ml Case of 36 24 mm

€25.20 As low
as
€20.32
plus VAT &
Shipping

Product Images



Description

Nalgene™ Unitary™ LDPE Syringe Bottles Dispense contents without tipping or shaking with Thermo Scientific™ Nalgene™ Unitary™ LDPE Syringe Bottles, which are faster and easier to fill with their wide neck. Capacity: 250 ml, 8 oz. Capacity (anglo): 8 oz. Color: Translucent Description: LDPE bottle/tube neck Quantity Per Carton: 36 Quantity Per Package: 4 Type: Rinse Bottle Closure Included: Yes Closure Size: 24 mm Closure Material: Polypropylene Closure Type: Screw Cap Graduated: No Height (Anglo-American): 5.42 in. Opening: Wide Shape: Cylindrical Dispensing through Closure: No Unit Size: Case of 36

Additional Information

No.	CD-0591
Manufacturer (Brand)	Nalgene
VGKL number	228724022
EAN	4058072165680
Gross weight	1.26
Transport temperature	Room temperature
Material	Low density polyethylene (LDPE)
Volume	250 ml
Country of origin	Mexico
Packaging volume	0.02700000 cc
Customs tariff number	39233010

