



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Unitary™ LDPE Squirt Bottles 500 ml Case of 24 38 mm

€30.64 As low
as
€24.68
plus VAT &
Shipping

Product Images



Description

Nalgene™ Unitary™ LDPE Syringe Bottles Dispense contents without tipping or shaking with Thermo Scientific™ Nalgene™ Unitary™ LDPE Syringe Bottles, which are faster and easier to fill with their wide neck. Capacity: 500 ml, 16 oz. Capacity (Anglo-American): 16 oz. Color: Translucent Description: LDPE bottle/tube neck Quantity Per Carton: 24 Quantity Per Package: 4 Type: Rinse Bottle Closure Included: Yes Closure Size: 38 mm Closure Material: Polypropylene Closure Type: Screw Cap Graduated: No Height (Anglo): 6.81 in. Opening: Wide Shape: Cylindrical Dispensing through Closure: No Unit Size: Case of 24

Additional Information

No.	CD-0592
Manufacturer (Brand)	Nalgene
VGKL number	228724023
EAN	4058072123161
Gross weight	1.98
Transport temperature	Room temperature
Material	Low density polyethylene (LDPE)
Volume	500 ml
Country of origin	Mexico
Packaging volume	0.03500000 ccm
Customs tariff number	39233010

