



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz automatic tipper complete, with dosing
attachment 5 ml, with 1,000 ml storage bottle,
NS 29/32**

**€53.08
plus VAT &
Shipping**

Product Images



Description

Automatic tipping machine

Dosing attachment with storage bottle (1 liter); ground joint: NS 29/32

Additional Information

No.	LZ-1448
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3451
EAN	4058072311872
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	5 ml
Material	Borosilicate glass
for bottles with volume	1000 ml
Spare part	LZ-1457
Material bottle	Borosilicate glass

