



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz pycnometer acc. to Gay-Lussac, not calibrated, nominal volume 25 cm<sup>3</sup>

**€33.98**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Gay-Lussac **pycnometer** not adjusted

Gay-Lussac pycnometer acc. to DIN ISO 3507 and supplementary to the standard, unadjusted, stopper NS 10/19 with capillary, upper stopper end ground and polished, nominal volume imprint on the bottom

## Additional Information

---

No.	LZ-1562
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3565
EAN	4058072313005
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass 3.3
Volume	25 ml
Form	pear-shaped
Quantity pistons	1
Standard ground sleeve	10/19
Standard ground core	10/19

