



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz pycnometer acc. to Gay-Lussac, not  
calibrated, nominal volume 100 cm<sup>3</sup>**

**€48.77  
plus VAT &  
Shipping**

## Product Images



## Description

---

Gay-Lussac **pycnometer** not adjusted

Gay-Lussac pycnometer acc. to DIN ISO 3507 and supplementary to the standard, unadjusted, stopper NS 10/19 with capillary, upper stopper end ground and polished, nominal volume imprint on the bottom

## Additional Information

---

No.	LZ-1564
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3567
EAN	4058072313029
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass 3.3
Volume	100 ml
Form	pear-shaped
Quantity pistons	1
Standard ground sleeve	10/19
Standard ground core	10/19

