



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® scintillation vials made of glass, 2 ml,  
58 x 27 mm screw cap cork-aluminum seal; PU:  
100 pcs./pack**

**€87.40  
plus VAT &  
Shipping**

## Product Images



## Description

---

The Ecolab scintillation vials are made of borosilicate glass 3.3.  
The scope of delivery includes screw caps with a cork-aluminum seal.

## Additional Information

---

No.	E-0065
Manufacturer (Brand)	neoLab
EAN	4058072082956
Transport temperature	Room temperature
Color	transparent (clear)
Volume	20 ml
Color bottle	transparent (clear)
Height	58 mm
Bottles properties	with closure/lid bottle autoclavable
Floor outside	flat
Floor inside	flat
DM bottle	27 mm
DM Neck	22 mm
Color closure/ lid	white
Thread	22-400 (US)
Neck size	far
Material bottle	Lime soda glass
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed

