




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Dropper bottle brown, screw pipette, 100 ml

€5.02 As low as

€4.77

**plus VAT &
Shipping**

Product Images



Description

Round dropper bottles with screw-pipette and rubber teat.

Made of amber glass for light-sensitive media.

Small quantities of liquid can be removed cleanly, dripped on and dispensed well.

Tightly closing.

Additional Information

No.	1-7244
Manufacturer (Brand)	neoLab
EAN	4058072046002
Color	brown
Volume	100 ml
Color bottle	brown
Type closure/ lid	Pipette dropper (screwable)
Bottles properties	with closure/lid bottle autoclavable
Material	Soda-lime glass (normal glass) 2
DM bottle	47.1 mm
DM Neck	18 mm
Color closure/ lid	white red
Neck size	narrow
Material bottle	Soda-lime glass (normal glass) 2
Material closure/ lid	Polypropylene (PP)

