



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® pressure pump sprayer, 1 L, with short nozzle

€24.90
plus VAT &
Shipping

Product Images



Description

For aqueous solutions, with unscrewable bottle.

Pump air into the container via a built-in pump lever (20 - 30 pump impacts depending on the filling level). With the spray button, the spray mist or jet can be dosed and easily dispensed. By adjusting the nozzle, it is possible to select between nebulization and spray jet. The spray intensity depends on the consistency of the medium and the overpressure.

Order No. 2-1985: with short nozzle

Color: black/white

Additional Information

No.	2-1985
Manufacturer (Brand)	neoLab
EAN	4058072059248
Transport temperature	Room temperature
Color	black white
Material	Polypropylene (PP) High-density polyethylene (HDPE) Polyamide (PA. non-specific) Steel
Volume	1000 ml
Gas type	Air

