



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Pressure pump atomizer 1 L, for solvents

€25.70
plus VAT &
Shipping

Product Images



Description

Order No. 2-1987 For isopropanol.

Pump air into the container via a built-in pump lever (20 - 30 pump pulses depending on filling level). With the spray button, the spray mist or jet can be dosed and easily dispensed. By adjusting the nozzle, it is possible to choose between nebulization and spray jet. The spray intensity depends on the consistency of the medium and the overpressure.

Additional Information

No.	2-1987
Manufacturer (Brand)	neoLab
EAN	4058072159061
Transport temperature	Room temperature
Color	black red white
Material	Polypropylene (PP) High-density polyethylene (HDPE) Polyamide (PA. non-specific) Steel Fluoroelastomer (Viton®)
Volume	1000 ml
Gas type	Air

