



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Pump spray bottle 10 ml, brown

**€2.85** As low as

**€2.72**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Lightweight vials with protective cap made of PS and riser tube made of HDPE.  
Ideal for spraying small quantities of liquid evenly and finely.

Pump sprayer with DIN 18 cap.

Spray angle 35 degrees, spray volume per pumping operation 0.05 ml.

## Additional Information

---

No.	1-7245
Manufacturer (Brand)	neolab
EAN	4058072046019
Color	brown
Material	Polystyrene (PS)
Volume	10 ml
Color bottle	brown
Type closure/ lid	Pump sprayer
Bottles properties	with closure/lid bottle autoclavable
Material	Polystyrene (PS)
Neck size	narrow
DM bottle	25 mm
DM Neck	18 mm
Color closure/ lid	transparent (clear)
Material bottle	Polystyrene (PS)
Material closure/ lid	Polystyrene (PS) Polypropylene (PP)

