




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Pump spray bottle 20 ml, brown

€3.09 As low as

€2.72

**plus VAT &
Shipping**

Product Images



Description

Lightweight vials with protective cap made of PS and riser tube made of HDPE.
Ideal for spraying small quantities of liquid evenly and finely.

Pump sprayer with DIN 18 cap.

Spray angle 35 degrees, spray volume per pumping operation 0.05 ml.

Additional Information

No.	1-7246
Manufacturer (Brand)	neolab
EAN	4058072046026
Color	brown
Material	Polystyrene (PS)
Volume	20 ml
Color bottle	brown
Type closure/ lid	Pump sprayer
Bottles properties	with closure/lid bottle autoclavable
Material	Polystyrene (PS)
Neck size	narrow
DM bottle	28.5 mm
DM Neck	18 mm
Color closure/ lid	transparent (clear)
Material bottle	Polystyrene (PS)
Material closure/ lid	Polystyrene (PS) Polypropylene (PP)

