



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Cryomaster® Cryotubes with external thread,  
2.0 ml, stand ring, 1D barcode, white screw cap,  
sterile, 20x25 pcs**

**€137.70  
plus VAT &  
Shipping**

## Product Images



## Description

---

The **Cryomaster® cryotubes with external thread** and stand ring are made of high-quality pure polypropylene. Thus, these cryovials are suitable for temperatures from -196°C to +121°C and can be used, among other things, for freezing samples with liquid nitrogen. A 1D barcode printed on the side of the tube makes these cryovials suitable for biobanking. The screw cap with external thread does not require an O-ring due to the precision fit, which further prevents potential contamination. The tubes are 20 x 25 sterile pouches and are DNase-, RNase-, endotoxin- and DNA-free.

## Additional Information

No.	4-6100
Manufacturer (Brand)	Cryomaster
EAN	4058072091279
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	12 mm
Color	White
Material	Polypropylene (PP)
Volume	2 ml
sterile	Yes
Height	49 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Cryotube Features	with closure/lid with labeling field graduated RNase-free RNA-free DNase-free DNA-free
Type of sterility	radiation sterile
Barcode	1D barcode lateral
Floor outside	flat. with stand ring
Floor inside	around
Type Tube/ Micro Vessel	Cryotube
suitable for cooling in	Nitrogen tank

