



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Cryotubes 1.2 ml, sterile, external
thread, stand ring, 100 pcs./pack**

**€50.50
plus VAT &
Shipping**

Product Images



Description

Cryotube (PP) with screw cap (PP) with silicone seal, 12.5 mm Ø.

With quick-release cap, autoclavable.

Temperature resistant to -196°C, only for storage in the vapor phase of liquid nitrogen. Graduated, with labeling field.

RNase-, DNase-, pyrogen-, DNA-free, gamma sterilized.

7-8051 is centrifugable up to 17000 g.

In resealable bags

Accessories: Lid inserts, colored (Order No. 7-8070 ff.)

Additional Information

No.	7-8050
Manufacturer (Brand)	neoLab
EAN	4058072042011
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	12,5 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	1.2 ml
sterile	Yes
Height	43 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Cryotube Features	pyrogen-free with closure/lid with seal with labeling field graduated color-codeable RNase-free DNase-free DNA-free
Floor outside	conical. with stand ring
Floor inside	around
Material gasket / septum	Silicone
suitable for cooling in	Nitrogen container Ultrat freezer Deep freezer Freezer / freezer
Type Tube/ Micro Vessel	Cryotube
centrifugable up to	17000 g

