



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene cryotubes 2.0 ml, 13.5 x 48.3 mm, 500 pcs/pack

€17.75 As low as
€14.25
**plus VAT &
Shipping**

Product Images



Description

Cryogenic tubes made of PP with screw cap, self-standing, for storage in the vapor phase of liquid nitrogen. Suitable for sample storage down to -196°C . White marking field, standardized size, non-cytotoxic, pyrogen-free, 13.5 mm \varnothing , radiation sterilized.
Packaging: 25 pcs. in bag, 20 bags/VE

Additional Information

No.	CD-1508
Manufacturer (Brand)	Nalgene
old neoLab article no.	3-6001
EAN	4058072178253
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	13,5 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	2 ml
sterile	Yes
Height	48.3 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Cryotube Features	pyrogen-free with closure/lid with labeling field graduated
Floor outside	conical. with stand ring
Floor inside	around
Type Tube/ Micro Vessel	Cryotube
suitable for cooling in	Nitrogen container Ultrat freezer Deep freezer Freezer/ freezer

