



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryovials, sterile, clear, 0.5 ml

€236.20
plus VAT &
Shipping

Product Images



Description

Cryovials, conical, sterile, clear, 0.5 ml, with lid.

- Material PP
- autoclavable at 121°C
- applicable down to -86°C
- usable in the vapor phase of liquid nitrogen
- screwed on lid
- gamma sterilized

Accessories for cryobox with hinged lid (order no. 2-0025,2-0030,2-0031)

Additional Information

No.	2-0020
Manufacturer (Brand)	neoLab
EAN	4058072047306
autoclavable at 121°C	Yes
diameter	7 mm
Color	transparent (clear)
Material	Polypropylene (PP)
Volume	0.5 ml
sterile	Yes
Type closure/ lid	Screw cap
Cryotube Features	with closure/lid with seal
Type of sterility	radiation sterile
Floor outside	conical
Floor inside	conical
Material gasket / septum	Silicone
TBST MAX	121 °C
TBST MIN	-86 °C
Type Tube/ Micro Vessel	Cryotube
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

