



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Cryotubes for Long Term Storage 1.5 ml Case of 500 Sterile

€27.00 As low
as
€22.00
plus VAT &
Shipping

Product Images

Description

Nalgene™ Cryotubes for Long-Term Storage Improve aseptic methods with Thermo Scientific™ Nalgene™ Cryotubes for Long-Term Storage with external threads, ideal for research applications as well as general laboratory storage. Single-handed operation: Yes Capacity (Anglo-American): 0.0507 oz. Closure material: Polypropylene Closure type: Screw cap DNase-free: Yes Packaging quantity: 500/carton Printed: Yes Graduated: Yes Type: Cryogenic screw cap Certification/Conformity: CE-marked for storage of diagnostic specimens in EU countries Description: Tubes, sterile, 1.5 ml Outer diameter (Anglo-American): 0.46 in. Autoclavable: Non-autoclavable Capacity (metric): 1.5 ml Thread Type: Outer Length (Anglo-American): 1.87 inches Quantity per case: 500 Quantity per pack: 25 RNase-free: Yes Product line: Nalgene™ Unit Size: Case of 500

Additional Information

No.	CD-1511
Manufacturer (Brand)	Nalgene
VGKL number	238751020
EAN	4058072219376
Gross weight	1.06
Transport temperature	Room temperature
Country of origin	USA
sterile	Yes
Packaging volume	0.010000 ccm
Customs tariff number	39269097

