



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Cryomaster® Cryotubes with internal thread,  
2.0 ml, stand ring, 1D barcode, white screw cap,  
sterile, 20x25 pcs**

**€137.70  
plus VAT &  
Shipping**

## Product Images



## Description

---

The **Cryomaster® cryotubes with internal thread** and stand ring are made of high-quality, ultra-pure polypropylene. Thus, these cryovials are suitable for temperatures from -196°C to +121°C and can be used, among other things, for freezing samples with liquid nitrogen. A 1D barcode printed on the side of the tube makes these cryovials suitable for biobanking. The screw cap with internal thread contains an O-ring to ensure maximum tightness.

The tubes are sterile pouched 20 x 25 and are DNase-, RNase-, endotoxin- and DNA-free.

## Additional Information

No.	4-6106
Manufacturer (Brand)	Cryomaster
EAN	4058072091996
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	11 mm
Color	White
Material	Polypropylene (PP)
Volume	2 ml
sterile	Yes
Height	49 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Cryotube Features	with closure/lid with internal thread with seal with labeling field graduated RNase-free RNA-free DNase-free DNA-free
Floor outside	flat. with stand ring
Floor inside	around
Material gasket / septum	Nitrile rubber (Buna-NO®)
Type of sterility	radiation sterile
Barcode	ID barcode lateral
suitable for cooling in	Nitrogen tank
Type Tube/ Micro Vessel	Cryotube

