



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® Cryotubes w. internal thread, sterile,  
self-standing, 4 ml, 100 pcs./pack**

**€58.50  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Cryotube (PP) with screw cap (PP) with silicone seal, 12.5 mm Ø.

With quick-release cap, autoclavable.

Temperature resistant to -196°C, only for storage in the vapor phase of liquid nitrogen. Graduated, with labeling field.

RNase-, DNase-, pyrogen-, DNA-free, gamma sterilized.

Round-bottomed tubes centrifugable up to 14000 g.

In resealable bags.

Accessories: Lid inserts, colored (Order No. 7-8070 ff.)

## Additional Information

No.	7-8064
Manufacturer (Brand)	neoLab
EAN	4058072042110
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	12,5 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	4 ml
sterile	Yes
Height	72 mm
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-196 °C
Cryotube Features	pyrogen-free with closure/lid with internal thread with seal with labeling field graduated color-codeable RNase-free DNase-free DNA-free
Floor outside	conical. with stand ring
Floor inside	around
Material gasket / septum	Silicone
suitable for cooling in	Nitrogen container Ultrat freezer Deep freezer Freezer/ freezer
Type Tube/ Micro Vessel	Cryotube
centrifugable up to	14000 g

