



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand disposable cuvettes PMMA, semi-micro,  
1.5-3.0 ml, 100 pcs./pack**

**€15.08** As low as  
**€12.56**  
**plus VAT &  
Shipping**

## Product Images

---



## Description

---

for kinetic measurements, nest numbered.

Cell: 12.5 x 12.5 x 45 mm

Window macro cell: 10 x 35 mm

Window semi-micro cell: 4.5 x 23 mm

Layer thickness: 10 mm

Polymethylmethacrylate (PMMA)

Application range from 300 nm

Standard deviations

at 320 nm  $\leq \pm 0.004$  Abs.

Packaging unit: 100 pcs./pack

## Additional Information

---

No.	BN-2095
Manufacturer (Brand)	Brand
old neoLab article no.	E-1643
EAN	4058072213220
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Volume MAX	3 ml
Volume MIN	1.5 ml
Wide	12.5 mm
Height	45 mm
Length	12.5 mm
Depth	12.5 mm
Cuvette properties	for 1xuse
Layer thickness	10.0 mm
Wavelength	from 300 nm
Standard deviation	at 320 nm $< \pm 0.004$
Type cuvette	Semi-micro

