



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® disposable cuvettes made of PS, semi-
micro, 1.5-3.0 ml, 100 pcs./pack**

**€9.53
plus VAT &
Shipping**

Product Images



Description

for kinetic measurements, nest-numbered.

Cell: 12.5 x 12.5 x 45 mm

Window macro cell: 10 x 35 mm

Window semi-micro cell: 4.5 x 23 mm

Layer thickness: 10 mm

Polystyrene (PS)

Application range from 340 nm

Standard deviations

at 360 nm $\leq \pm 0.005$ Abs.

Packaging unit: 100 pcs./pack

Additional Information

No.	E-1642
Manufacturer (Brand)	neoLab
EAN	4058072066178
Transport temperature	Room temperature
Material	Polystyrene (PS)
Volume MAX	3 ml
Volume MIN	1.5 ml
Wide	12.5 mm
Height	45 mm
Length	12.5 mm
Depth	12.5 mm
Cuvette properties	for 1xuse
Layer thickness	10.0 mm
Wavelength	from 340 nm
Standard deviation	at 300 nm $< \pm 0.005$
Type cuvette	Semi-micro

