



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma Half Micro Cell 117.004F-QS, SD 10x4mm, VOL 1400µl

€514.00
plus VAT &
Shipping

Product Images



Description

Sealable half-micro cell type 117.104F-QS with ISO thread GL 14 and open screw cap with silicone seal (septum), for fluorescence measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 1500 µl, path length 10 x 4 mm, outer dimensions 56 x 12.5 x 12.5 mm, inner width 4 mm, bottom thickness 1.25 mm

Additional Information

No.	KI-0562
Reg. trademark	SUPRASIL®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	105511174
EAN	4058072173258
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01280 kg
Material	Quartz glass
Volume MAX	1.4 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.06900 m
Packing height	0.02200 m
Packing depth	0.09400 m
Customs tariff number	90279050
Cuvette properties	for fluorometry
Wavelength	200 to 2500 nm

