



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand cuvette, 4 optical windows, macro, nest
number sorted, UV-transparent, minimum
filling volume 2.5 ml**

**€41.80
plus VAT &
Shipping**

Product Images



Description

Disposable cuvettes 4 optical windows, macro, PS or UV transparent, nest number sorted, minimum filling volume 2.5 ml. Suitable for fluorescence spectroscopy. The UV version can be used from a wavelength of 230 nm and shows no autofluorescence. Quality features:

Optically flawless transmission range Manufacture under controlled room conditions and fully automated packaging Sorted according to nest numbers to ensure the lowest possible deviations in absorbance values The UV-transparent cuvettes replace sensitive and expensive glass or quartz glass cuvettes in many areas

Additional Information

No.	BN-2096
Manufacturer (Brand)	Brand
VGKL number	104975125
EAN	4058072213237
HSNumber	90279000
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.23100 kg
Country of origin	Germany
Customs tariff number	90279000

