



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hellma macro cuvette 100-QS SD 10mm, VOL 3500µl

**€96.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Macro Cell Type 100-QS with folded lid for absorbance measurements, window material High Performance Quartz Glass (200 nm - 2500 nm), volume 3500 µl, path length 10 mm, outer dimensions 45 x 12.5 x 12.5 mm, inner width 9.5 mm, bottom thickness 1.5 mm

## Additional Information

---

No.	KI-0432
Reg. trademark	SUPRASIL®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	195511102
EAN	4058072176501
Gross weight	0.02
Transport temperature	Room temperature
Net weight	0.00740 kg
Material	Quartz glass
Volume MAX	3.5 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	10.0 mm
Wavelength	200 to 2500 nm

