



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma micro cuvette 104.002-OS SD 10mm, VOL 700µl

€84.60
plus VAT &
Shipping

Product Images



Description

Micro cuvette type 104.002-OS with folded lid for absorbance measurements, window material optical special glass (320 nm - 2500 nm), volume 700 µl, path length 10 mm, outer dimensions 45 x 12.5 x 12.5 mm, inner width 2 mm, bottom thickness 3.2 mm

Additional Information

No.	KI-0466
Reg. trademark	OS ®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	195514161
EAN	4058072181666
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01660 kg
Material	Optical glass
Volume MAX	0.7 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	10.0 mm
Wavelength	320 to 2500 nm

