



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hellma micro cuvette 104.002F-QS SD 10x2mm, VOL 700µl

**€346.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Micro cuvette type 104.002F-QS with folded lid for fluorescence measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 700 µl, path length 10 mm x 2 mm, outer dimensions 45 x 12.5 x 12.5 mm, inner width 2 mm, bottom thickness 1.25 mm

## Additional Information

---

No.	KI-0472
Reg. trademark	SUPRASIL®
Manufacturer (Brand)	Hellma
VGKL number	195514182
Dual use indicator	false
EAN	4058072174651
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01390 kg
Material	Quartz glass
Volume MAX	0.7 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for fluorometry
Wavelength	320 to 2500 nm

