



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hellma micro cuvette 108.002-QS SD 10mm, VOL 500µl

**€191.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Micro cuvette type 108.002-QS with folded lid for absorbance measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 500 µl, path length 10 mm, outer dimensions 45 x 12.5 x 12.5 mm, inner width 2 mm, bottom thickness 9 mm

## Additional Information

---

No.	KI-0512
Reg. trademark	SUPRASIL®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	105514192
EAN	4058072785673
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01420 kg
Material	Quartz glass
Volume MAX	0.5 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	10.0 mm
Wavelength	200 to 2500 nm

