



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hellma micro cuvette 115-QS SD 10mm, VOL 400µl

**€198.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Micro cuvette type 115-QS with stopper made of PTFE for absorbance measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 400 µl, path length 10 mm, outer dimensions 40 x 12.5 x 12.5 mm, inner width 2 mm, bottom thickness 1.25 mm

## Additional Information

---

No.	KI-0555
Reg. trademark	SUPRASIL®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	195514192
EAN	4058072785963
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01260 kg
Material	Quartz glass
Volume MAX	0.4 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	10.0 mm
Wavelength	200 to 2500 nm

