



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma micro cuvette 115B-QS SD 10mm, VOL 400µl

€373.00
plus VAT &
Shipping

Product Images



Description

Micro cuvette type 115B-QS with stopper made of PTFE for absorbance measurements, window material High Performance Quartz glass (200 nm - 2500 nm), side walls and bottom black, volume 400 µl, path length 10 mm, outer dimensions 40 x 12.5 x 12.5 mm, inner width 2 mm, bottom thickness 1.25 mm

Additional Information

No.	KI-0556
Reg. trademark	SUPRASIL®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	105514195
EAN	4058072785970
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01290 kg
Material	Quartz glass
Volume MAX	0.4 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	10.0 mm
Wavelength	200 to 2500 nm

