



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma micro cuvette 115F-QS SD 10x2mm, VOL 400µl

€372.00
plus VAT &
Shipping

Product Images



Description

Micro cuvette type 115F-QS with stopper made of PTFE for fluorescence measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 400 µl, layer thickness 10 x 2 mm, outer dimensions 40 x 12.5 x 12.5 mm, inner width 2 mm, bottom thickness 1.25 mm

Additional Information

No.	KI-0557
Reg. trademark	SUPRASIL®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	105514184
EAN	4058072785987
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01300 kg
Material	Quartz glass
Volume MAX	0.4 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for fluorometry
Wavelength	200 to 2500 nm

