



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf UVette® Starter Set, 50 – 2,000 µL

€102.00
plus VAT &
Shipping

Product Images



Description

UVette® Starter Set, 80 UVettes and 1 universal adapter for light path height 15 mm, convertible to 8.5 mm, 50 – 2,000 µL – Individually packed – Quality PCR clean and protein-free – UV-transparent plastic, free of fluoropolymers and other halogenated hydrocarbons – Two optical path lengths in a single cuvette: To switch between 10 mm and 2 mm, simply rotate the UVette by 90° – Volume markings at 500 µL and 1.000 µL – Individually packed and according to certified quality PCR clean and protein-free for sensitive (e.g. RNA) and precious samples – Deep optical window prevents scratches – Funnel-shaped cuvette bottom for optimal filling, roughened grip surface for labeling – Ideally suited for use in the Eppendorf BioPhotometer and Eppendorf BioSpectrometer, adapters are available for use in all other conventional spectrophotometers. – Routine pack for convenient access to each UVette and for safe storage in a resealable box

Additional Information

| | |
|-----------------------|--------------------|
| No. | VB-0986 |
| Manufacturer (Brand) | Eppendorf |
| VGKL number | 296991007 |
| EAN | 4058072411268 |
| Gross weight | 0.316 |
| HSNumber | 90279000 |
| Transport temperature | Room temperature |
| Storage temperature | Room temperature |
| Net weight | 0.28500 kg |
| Material | Plastic |
| Country of origin | Germany |
| Packing width | 0.31000 m |
| Packing height | 0.10000 m |
| Packing depth | 0.13000 m |
| Packaging volume | 5000.000000 ccm |
| Customs tariff number | 90279000 |
| Wavelength | 220 nm to 1,600 nm |

