



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand UV cuvette micro, center height 15 mm,
filling volume 70 µl to 550 µl**

€132.75 As low
as
€119.75
plus VAT &
Shipping

Product Images



Description

Specially developed for the photometric determination of proteins, ssDNA, dsDNA, RNA and oligonucleotides in the UV range. Applicable for wavelengths from 230 nm, ideally suited for measurements at 260 nm, 280 nm and in the visible wavelength range. Standard layer thickness 10 mm. Even small sample quantities from 70 µl are sufficient for measurement. Individually packaged UV cuvettes micro are DNase, DNA and RNase free

Additional Information

No.	BN-2104
Manufacturer (Brand)	Brand
VGKL number	294959230
EAN	4058072213312
Gross weight	1.2019
HSNumber	90279000
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.73450 kg
Country of origin	Germany
Packing width	0.22700 m
Packing height	0.20500 m
Packing depth	0.36000 m
Packaging volume	16752.60000000 ccm
Customs tariff number	90279000

